

DISPOSABLE INFLATION DEVICE : AGI-40X




Features




- Simple design
- Maximum allowable pressure : 40 atm
- One touch release lever
- Maximum liquid volume : 20 ml



Catalog No.	Product name	Description	Configuration
-------------	--------------	-------------	---------------

AGI-40X	Disposable Inflation Device	To inflate/deflate balloon catheters as well as to monitor pressure within the balloon	
---------	-----------------------------	--	--



- Read instruction prior to use

Marketing Approval number : 13B3X10277110320

Manufactured by : AG and International Inc.

3-20-4, Hongo, Bunkyo-ku, Tokyo 113-0033 Japan

Tel : +81-3-6801-6101 Fax : +81-3-6273-0170

www.ag-intl.co.jp

Indication

Vascular case

- Cardiovascular intervention case
- Percutaneous Transluminal Angioplasty
- Aortic valve stenosis case (TAVI)
- Vascular Access Interventional Therapy
- Endovascular Aneurysm Repair
- Thoracic Endovascular Aneurysm Repair

Non-vascular case (Interventional radiology)

- B-RTO: balloon occluded retrograde transvenous obliteration
- PTO: Percutaneous Transhepatic Obliteration
- PVP: Percutaneous Vertebroplasty)
- BKP: Balloon Kyphoplasty
- VBS: Vertebral Body Stenting

conclusion

AGI-40X inflation device allows 40atm maximum allowable pressure and it increases to treat more lesions, not only vascular stenotic lesion, but also non-vascular atenotic lesions