

Performing by design

CTO balloon "SHOW" 翔

PTCA Balloon Catheter

Designed to cross the tightest lesions



Committed to excellence in developing Arterial Therapeutic Systems

CTO balloon "SHOW" 翔

PTCA Balloon Catheter

- > Unsurpassed crossability and pushability.
- > Lowest crossing profile.
- > Innovative shaft material for OTW version.


		Diameter (mm)			
		1.00	1.10	1.25	1.50
Length (mm)	8	Available	Available	Available	Available
	10	Available	Available	Available	Available
	15	Available	Available	Available	Available
	20	Available	Available	Available	Available
	25	Available	Available	Available	Available

		Pressure (bars)									
		2	4	6*	8	10	12	14	16	18**	20
Balloon diameter (mm)	ø1.00	0.96	0.98	1.00	1.02	1.04	1.06	1.08	1.10	1.12	1.15
	ø1.10	1.06	1.08	1.10	1.12	1.14	1.17	1.19	1.21	1.23	1.27
	ø1.25	1.20	1.23	1.25	1.28	1.30	1.32	1.35	1.38	1.40	1.44
	ø1.50	1.44	1.47	1.50	1.53	1.56	1.59	1.62	1.65	1.68	1.73


Available

* Nominal pressure ** Rated Burst Pressure

Sirius CTO RX



Sirius CTO OTW



- > Guidewire compatibility: 0.014" maximum recommended
- > Guiding catheter compatibility: 5F (0.058" ID)



AAA INTL Ltd

33 Bridge Street, Hereford, HR4 9DQ UK

Tel +44 1926 800153

E-mail : info@aaa-intl.com

Web-site : http://www.aaa-intl.com